

WD-TEC

Fuego<sup>®</sup> SCS



Safety Laboratory  
Gas Burner



DOMINIQUE FUSCHER

# Innovative technology that thinks with you!

The new Fuego SCS - Series.

**Graphical installation instructions**

**Zero-pressure shut-off**

**Gas consumption display**

**Acoustic signals as operating aids**

**User-friendly stainless steel knobs**

**Animated, high-contrast graphic display**

**Individually adjustable Single or DoubleClick IR motion sensor**

**Automatic unit cut-off**

**Cooling reminder for inoculation loops**

**Temperature regulation for heating media**

DOMINIQUE

## Fuego<sup>®</sup> pro

### Battery performance

The exclusive Fuego SCS pro provides maximum flexibility. Wireless operation using state-of-the-art rechargeable battery technology guarantees independent operation. Two standard rechargeable batteries allow up to 9 hours of continuous operation. This corresponds to approximately 2000 inoculation loop flaming operations. The Fuego SCS pro has an integrated fast charging function and can be recharged in just 3 hours. The range of functions is identical to those of the Fuego SCS.



	Fuego <sup>®</sup> basic	Fuego <sup>®</sup> basic RF	Fuego <sup>®</sup>	Fuego <sup>®</sup> pro
<b>Programs</b>				
Button function	✓	✓	✓	✓
Pedal Standard	✓	✓	✓	✓
Pedal Start-Stop	✓	✓	✓	✓
Sensor Auto-Off	✓	✓	✓	✓
Sensor Start-Stop	✓	✓	✓	✓
<b>NEW: Temperature regulation</b>				
<b>Safety Control System</b>	✓	✓	✓	✓
with acoustic signals	✓	✓	✓	✓
Residual heat display	✓	✓	✓	✓
<b>NEW: Animated graphic display</b>				✓
Foot pedal included	✓	✓ (radio)		
Single- or DoubleClick IR-Sensor			✓	✓
Battery operation with tilt sensor				✓
<b>NEW: Cooling reminder</b>			✓	✓
<b>NEW: Zero-pressure shut-down</b>			✓	✓
<b>NEW: Gas consumption display</b>			✓	✓
<b>NEW: 2 user accounts</b>			✓	✓
Gas inlet filter	✓	✓	✓	✓
Removable & decomposable burner head	✓	✓	✓	✓
Housing, burner head, gas connection & controls made of stainless steel	✓	✓	✓	✓
DVGW approved	✓	✓	✓	✓
<b>Art.-No.:</b>	<b>8.201.000</b>	<b>8.203.000-RF</b>	<b>8.200.000</b>	<b>8.204.000</b>

Removable & decomposable burner head, dishwasher-safe, autoclavable, made completely of stainless steel and ceramics



Burner shaft can be opened at the flick of a wrist for cleaning purposes



Tilt adjustment



Holding device for inoculation loop holders

DOMINIQUE

SCHER