

FlyStuff Flowbuddy™ Flow Regulator

The Flowbuddy Flow Regulator adds convenience and precision flow control to any CO₂ Drosophila workstation.



Wall-Mount Flowbuddy

Benchtop Flowbuddy, Complete

Cat #	Description
59-122B	Benchtop Flowbuddy, Alone
59-122BC	Benchtop Flowbuddy, Complete*
59-122BCU	Benchtop Flowbuddy, Complete* w/ Ultimate Flypad
59-122W	Wall-Mount Flowbuddy, Alone
59-122WC	Wall-Mount Flowbuddy, Complete*
59-122WCU	Wall-Mount Flowbuddy, Complete* w/ Ultimate Flypad

* Complete Flowbuddy system includes Flowbuddy flow regulator; blowgun with coiled tube and 0.05 x 1 3/16" needle (0.72 x 1 3/4" can be substituted free upon request); standard size flypad (large size can be substituted free upon request); and 6 feet 1/8" I.D. tubing.

FlyStuff Foot Valve

The FlyStuff Foot Valve is a simple and efficient tool designed to conserve CO₂ and reduce the risk of over-anesthetizing flies. When compared to a standard continuous-flow setup, the addition of this simple foot valve to your Drosophila workstation can conserve CO₂ by 500% or more. This can translate into huge savings. In addition to saving on CO₂ cost, the FlyStuff Foot Valve saves the time and tedium involved in having to frequently change out CO₂ tanks. The solid brass barbed input and output fit 1/8" I.D. tubing.



Foot Valve Complete System



Input /Output

Cat #	Description
59-121	Flystuff Foot Valve
59-121C	Flystuff Foot Valve, Complete*

* The FlyStuff Complete foot valve system includes: 1 foot valve; 1 Flypad; 1 Flypad Frame; 5 Feet 1/8" I.D. Tubing; 1 "T" Fitting; 1 blowgun with coiled tube and 0.05 x 1 3/16" needle (0.72 x 1 3/4" can be substituted free upon request.)



Luer lock for needles

Cat # 54-119N (Actual Size)

Cat # 54-119 (Actual Size)

Cat # 54-119L (Actual Size)

Flystuff Mini Blowgun w/ Cat # 54-119L Long Replacement Needle attached

FlyStuff Blowgun

The FlyStuff pistol grip blowgun offers variable flow control and is lightweight, ergonomic and the perfect complement to any CO₂ anesthetizing apparatus. The output is a 0.052" stainless steel needle tip with a nickel plated brass base. The needle length measures 1 3/16" in length (ideal for use with cotton, rayon or Flugs). The barbed input is a solid brass 1/8" I.D. and fits our standard size tubing.

Cat #	Description
54-104	Flystuff Blowgun
54-104M	Flystuff Blowgun Mini
54-119	Replacement Needle, 0.05 x 1 & 3/16"
54-119L	Replacement Needle, 0.05 x 1 & 3/4"
54-119N	Luer Lock Adapter

>500% CO₂ Savings!

