

Certificate of analysis nº: 24/00017765

Emission date	2024-10-14		
Reference	920085200813		
Description	VENABIO FLOW 60shA (07,90)x(12,70)		
Length	50ft		
Hose Batch	IS00125859	Material Batch	IS00125643
Manufacturing date	27.09.2024	Expiration date (According ISO 2230 "Guidelines for storage")	25.09.2034

Properties	Unit	Inspection Method	Limits	Results
Hardness	Shore A	ISO 868	60 +/- 5	59.0
Tensile strength	N/mm2	ISO 37	>7	10.2
Elongation	%	ISO 37	>300	640.0
Tear strength B	kN/m	ISO 34-1	>20	55.0
Density	g/cm3	ISO 2781	1,16 +/- 0,02	1.16

Certified Information

Process	Quality management system certified ISO 9001:2015, ISO 14001:2015 and EMAS. Produced in ISO7 classified clean room according to ISO14644.
---------	--

Product	FDA Standard 21 CFR 177.2600
	BfR Standard part XV
	USP In Vivo <88> Class VI at 121°C
	ISO 10993-5 & USP In vitro <87>
	ISO 10993-4, 6 & 10
	European Pharmacopoeia 3.1.9

Approval

Quality Department

Venair Ibérica S.A.U.
C/Cerdanya 26, Pol. Ind. Nord - 08226 - Terrassa (Barcelona) - Spain

It is the user's responsibility to ensure the suitability and safety of the VENAIR products for all intended uses. All the tests must be conducted in accordance with applicable regulatory requirements to determine the safety and effectiveness for use of the hoses in any particular application. User assumes all other risk, if any, including the risk of injury, loss or damage, direct or consequential, arising out of the use, misuse, or inability to use, this product. Vena, Venair and the Venair logo are trademarks of Venair Ibérica SA.