

## Certificate of analysis n°: 26/00005244

Emission date	2026-02-23		
Reference	760321000019		
Description	TECHNOSIL USP VI D3/4 (19,1 mm)		
Length	20m		
Hose Batch	123453	Material Batch	
Manufacturing date	19.09.2023	Expiration date (According ISO 2230 "Guidelines for storage")	

Properties	Unit	Inspection Method	Limits	Results
Hardness	Shore A	ISO 868	70 +/- 5	
Tensile strength	N/mm2	ISO 37	≥5	
Elongation	%	ISO 37	>300	
Tear strength B	kN/m	ISO 34-1	>30	
Density	g/cm3	ISO 2781	1,19 +/- 0,03	

## Certified Information

Process	Quality management system certified ISO 9001:2015, ISO 14001:2015 and ISO 13485:2016. Produced in ISO7 classified clean room according to ISO14644 in Venair's Spanish manufacturing plant.
---------	--

Product	FDA Standard 21 CFR 177.2600 Section (e), (f)
	BfR Standard part XV
	USP In Vivo <88> Class VI at 121°C
	USP In vitro <87>
	ISO 10993-4, 5, 10, 11 & 23
	Regulation (EU) No. 10/2011

## Approval

Quality Department

Venair Biotech S.L.U.

C/ l'Alguer 28, Pol. Ind. Nord - 08226 - Terrassa (Barcelona) - Spain

It is the user's responsibility to ensure the suitability and safety of the VENAIR products for all intended uses. All the tests must be conducted in accordance with applicable regulatory requirements to determine the safety and effectiveness for use of the hoses in any particular application. User assumes all other risk, if any, including the risk of injury, loss or damage, direct or consequential, arising out of the use, misuse, or inability to use, this product. Vena, Venair and the Venair logo are trademarks of Venair Ibérica SA.