

TouchNTuff® DermaShield® 73-701

PRODUCT PROTECTION

NE NEOPRENE

SPLASH

| STYLE # | LINER MATERIAL | COATING MATERIAL | GRIP DESIGN | CUFF STYLE | THICKNESS MM | POWDER-FREE | AQL | COLOUR | EN SIZE | LENGTH MM | PACKAGING |
|---------|----------------|------------------|-------------------|------------|--------------|-------------|-----|--------|---------------------------|-----------|-----------|
| 73-701 | N/A | Neoprene | Bisque fingertips | Straight | 0.175 | Yes | 1.0 | Green | 6, 6.5, 7, 7.5, 8, 8.5, 9 | 310 | See below |



DESCRIPTION

- Made of poly-chloroprene, TouchNTuff® is the most technically advanced sterile glove available today.
- Free of latex proteins and accelerators, it minimises both Type I and Type IV allergies associated with natural rubber.
- While TouchNTuff® has the comfort and flexibility of natural rubber gloves, at a similarly competitive price, its neoprene composition offers better chemical resistance. TouchNTuff®: for optimal product and worker protection.
- The TouchNTuff® 73-701 is sterile.

REMARKS

- Complies to EN455 parts 1 and 2



CATEGORY III

The most advanced glove solution in allergy prevention

PRIMARY INDUSTRIES



IDEAL APPLICATIONS

- Aseptic drug manufacturing
- Biotechnologies
- Medical device manufacturing
- Handling of cytostatics

PACKAGING

Pair pack: The glove is packaged one pair per poly wallet, sealed in a poly envelope; 10 pairs per sealed polybag; 5 sealed polybags per master bag; 4 master bags per case liner; 200 pairs per carton. All packaging is printed with cleanroom-compatible IPA-resistant ink. The polybags have an easy-tear opening. The gloves are triple-bagged for superior contamination control.

Print date: 09-04-2017

Ansell Healthcare Europe N.V. (European Head Office)

Riverside Business Park, Block J
 Boulevard International 55, 1070 Brussels, Belgium
 Tel. +32 (0) 2 528 74 00 • Fax +32 (0) 2 528 74 01 • Fax Customer Service +32 (0) 2 528 74 03
<http://www.ansell.eu> • E-mail info@ansell.eu