



23 ml	AR	III	Clear	-	611.1775145105.9F005
Quantity	Material	Hydr.cl.	Colour	Remark	Article number

Demands/Remarks :

Standard Tolerance	Quality	Dimensions				Eccentricity		Radius	Colour Break point :
		0-1	>1-60	>60-120	>120	<14.75	P 14.75	min: 0,0	Printing paint :
	fine	0,03	0,30	0,50	0,80	0,60	0,80	max: 3,0	Stress relieving
	middle	0,05	0,50	0,80	1,20	-1,00	-1,20		by temp off : 560 °C
All dimensions, specification and symbols according ISO 9187 - 1.		Changes		Name	Date	Drawn	Name	Date	Weight :
Scale : x:1				pk	151208	Control	hge	120103	xx,x ± 0,3 gr
							ad/hge	120103	Client : Stock

Name : Screw neck tubes, 23 ml., thread capalu 16, clear sodalime glass 9.1, type 3, dim.  $\varnothing 17.75 \times 144.5 \times 1.05$  mm., unprinted, according to NAFVSM 611.1775145105.9F005.01-cl-3.

Drawing no. : 611.1775145105.9F005.01

**Plant Nijmegen**  
 Marialaan 108, NL 6541 RP Nijmegen, Netherlands  
 Tel +31-24-3774670  
 www.naf-vsm.com, e-mail: info.nl@naf-vsm.com

**Plant Andeville**  
 Zone Industrielle, F 60570 Andeville, France  
 Tel +33-3-44.52.64.90  
 www.naf-vsm.com, e-mail: info.fr@naf-vsm.com



h:\sketch2\detail\00000000.000.001.00.130101