



Frickenhausen, 27<sup>th</sup> February 2019

## Product Information: Decapper-suitable Cryo.s™ with external thread and Cryo.s™ 48-way CapCarrier

Dear Sir or Madam,

We are continuously striving to optimise our product portfolio and to expand our product portfolio with attractive products.

Therefore we would like to inform you about the following products:

### Decapper-suitable Cryo.s™ with external thread

As of now, Decapper-suitable cryo tubes (2 ml and 4 ml) with external thread are available. The new cryo tubes will be available as new versions / variants of 126263-2DG and 127263-2DG.

In addition to the well-established tube design with starfoot base, external thread and barcoding, decapper-suitable Cryo.s™ feature an **internal contour in their screw cap**. This feature allows for a **precise interaction of automated decapping systems with the tube caps and an exact torque value to be applied for use within machine-assisted, high-throughput workflows**.

Please note that Cryo.s™ with internal thread (122263-2DG and 123263-2DG) have been equipped with a Decapper-suitable screw cap for several years now and can therefore also be used with decappers.

### Cryo.s™ 48-way CapCarrier

Furthermore, we are pleased to inform you that we will offer a **48-way CapCarrier as of now**.

The novel Cryo.s™ 48-way CapCarrier serves as a perfect receptacle for the **intermediate storage of screw caps** (e.g. after decapping, during tube filling, before recapping), providing a depot for 48 externally or internally threaded Cryo.s™ screw caps.

The following items are available **as of now**:

Order no.	Product Description
126263-2DC	CRYO.S, 2 ML, PP, ROUND BOTTOM, EXTERNAL THREAD, NATURAL SCREW CAP FOR DECAP., WRITING AREA, STARFOOT, NATURAL, STERILE, 100 PCS./BAG, DATAMATRIX ECC 200, 14X14 DOTS, BARCODE 128
127263-2DC	CRYO.S, 4 ML, PP, ROUND BOTTOM, EXTERNAL THREAD, NATURAL SCREW CAP FOR DECAP., WRITING AREA, STARFOOT, NATURAL, STERILE, 50 PCS./BAG, DATAMATRIX ECC200, 14X14 DOTS, BARCODE 128
803507	CAP CARRIER, 48 WAY, PC, BLACK, FOR STORAGE OF CRYO.S SCREW CAPS, STERILE, 1 PC./BAG

For detailed product information please refer to the attached document F074072.

If you are interested in these products, please contact your Export Manager, local subsidiary or GBO distributor.

Kind regards,

i.V. Michael Wirth  
Director Marketing & Product Management