

Autoclavable Bags BIOHAZARD

Autoclavable bags BIOHAZARD – Warning Instruction and sterilization indicator patch

Polypropylene (PP).

Autoclavable disposable waste bag of very strong PP-film for the sterilization and disposal of hazardous and contaminated labware.

- rectangle white indicator patch field turns into dark when autoclaved
- with BIOHAZARD symbol and safety imprint in English, German, French and Spanish language
- autoclavable at 134 °C, 20 min, for autoclaving, the bag must not be completely closed



W x H mm	Volume approx. l	Film Thickness	Packaging	Quantity per Pack	Order No.
300 x 500	6	50 µm	cardboard box 1 x 100	100	70 02 105
300 x 500	6	50 µm	cardboard box 1 x 500	500	70 02 100
400 x 780	22	50 µm	cardboard box 1 x 100	100	70 02 155
400 x 780	22	50 µm	cardboard box 1 x 500	500	70 02 150
600 x 800	32	50 µm	cardboard box 1 x 100	100	70 02 205
600 x 800	32	50 µm	cardboard box 1 x 500	500	70 02 200
700 x 1100	72	50 µm	cardboard box 1 x 75	75	70 02 305
700 x 1100	72	50 µm	cardboard box 1 x 350	350	70 02 300

Autoclavable bags BIOHAZARD red – Warning Instruction and sterilization indicator patch

Polypropylene (PP).

Distinctive red autoclavable bags for easy identification, for the disposal of contaminated labware.

- rectangle white indicator patch field turns into dark when autoclaved
- with BIOHAZARD symbol and safety imprint in German, English and French language
- autoclavable at 134 °C, 20 min, for autoclaving, the bag must not be completely closed



W x H mm	Volume approx. l	Film Thickness	Packaging	Quantity per Pack	Order No.
300 x 500	6	50 µm	cardboard box 1 x 500	500	70 03 200
400 x 780	22	50 µm	cardboard box 1 x 500	500	70 03 210
600 x 800	32	50 µm	cardboard box 1 x 500	500	70 03 220
700 x 1100	72	50 µm	cardboard box 1 x 350	350	70 03 230

Autoclavable bags BIOHAZARD yellow – Warning Instruction and sterilization indicator patch

Polyethylene (HD-PE).

Autoclavable disposable waste bag of very strong PE-film for the sterilization and disposal of hazardous and contaminated labware.

- rectangle white indicator patch field turns into dark when autoclaved
- with BIOHAZARD symbol and safety imprint in German, English and French language
- autoclavable at 121 °C, 20 min, for autoclaving, the bag must not be completely closed



W x H mm	Volume approx. l	Film Thickness	Packaging	Quantity per Pack	Order No.
300 x 500	6	50 µm	cardboard box 1 x 500	500	70 03 100
400 x 780	22	50 µm	cardboard box 1 x 500	500	70 03 110
600 x 800	32	50 µm	cardboard box 1 x 500	500	70 03 120
700 x 1100	72	50 µm	cardboard box 1 x 350	350	70 03 130