


# VIRUS COLLECTION AND TRANSPORT KITS

DESIGNED FOR VIRAL SPECIMEN COLLECTION, TRANSPORT AND PRESERVATION, SUCH AS FLU, BIRD FLU AND HFMD, ETC.



- Transport tube : 5ml / 10ml
- Transport medium : 1ml/tube , 3ml/ tube
- Sterile swab : Nylon applicator



**CITOTEST®**



1. Peel off the pouch and collection
2. Put the swab into the medium
3. Snap the stick of swab
4. Screw the tube cap tightly
5. Label the specimen information
6. Put into the Biohazard bag

- Disposable collection swab, composed of snappable plastic stick and nylon applicator .
- VTM Transport media can be chosen as needed.
- Ready to use and easy to tear package, effectively avoid cross contamination.
- Supplied with Biohazard specimen bag, ensure transport safe and reliable.

PACKING: 50set/box, 400sets/carton

P/N

DESCRIPTION

P/N	DESCRIPTION
2118-1501-99	Virus collection set, 5ml tube with 1ml VTM medium, one nasal nylon swab, one biohazard specimen bag
2118-1502-99	Virus collection set, 5ml tube with 3ml VTM medium, one nasal nylon swab, one biohazard specimen bag
2118-1503-99	Virus collection set, 10ml tube with 1ml VTM medium, one nasal nylon swab, one biohazard specimen bag
2118-1504-99	Virus collection set, 10ml tube with 3ml VTM medium, one nasal nylon swab, one biohazard specimen bag
2118-1552-99	Virus collection set, 5ml tube with 1ml VTM medium, one oral nylon swab, one biohazard specimen bag
2118-1517-99	Virus collection set, 5ml tube with 3ml VTM medium, one oral nylon swab, one biohazard specimen bag
2118-1553-99	Virus collection set, 10ml tube with 1ml VTM medium, one oral nylon swab, one biohazard specimen bag
2118-1513-99	Virus collection set, 10ml tube with 3ml VTM medium, one oral nylon swab, one biohazard specimen bag
2118-1514-99	Virus collection set, 10ml tube with 3ml VTM medium, one nasal nylon swab, one oral nylon swab, one biohazard specimen bag



**CITOTEST LABWARE MANUFACTURING CO., LTD**

Address: No.339, Beihai West Road, Haimen 226100, Jiangsu Province

Contact: TEL +86 513 82259348 ; FAX+86 513 82110195

Web: [www.citotest.com](http://www.citotest.com)

**CITOTEST®**