

# ONE-STEP

## GL45 Filtration Adaptor



### ✓ Patented 3-way connection design (M 429664)

Unique filtration adaptor can be connected to vacuum source, filter holder and GL45 storage bottle. (e.g., 1000ml, 2000ml, 5000ml) at the same time.

### ✓ Compatible with various kinds of filter holders (with No.8 silicone stopper)

Filtration adaptor can be connected with glass, stainless steel and PES filter holders with No.8 silicone stopper for a wide range of applications.

### ✓ Preventing cross-contamination by one-step filtration

Ideal for keeping the filtrated solution directly in a GL45 storage bottle without the need to be transferred from the filtering flask.

### ✓ Resistant to organic solvent, autoclavable

Filtration adaptor is made of PP material with PTFE O-ring having excellent corrosion resistance to a wide range of chemicals.

#### Ordering Information

197000-65 Filtration Adaptor

