

BioDolphin Suction Kit



Video Introduction



BioDolphin System

Introduction

BioDolphin's ergonomic handle and friendly designed lock switch can keep continuous suction without pressing button all the time. The innovative design of 8-channel adapter enables tip ejection with single hand, more convenient to use. The complete kit except rack are available for autoclave.

Applications

- Cell culture
- Waste suction
- ELISA wash solution

Benefits:

- User-friendly lock switch design
- Single-handed tip ejection
- Multi-functional storage rack
- Quick fitting change



The following color options are available for order: Light blue, green, pink and deep blue.

Complete kit

- 1 Handle
- 2 1-channel tip adaptor
- 3 1-channel S.S adaptor (80mm)
- 4 1-channel S.S adaptor (40mm)
- 5 8-channel S.S adaptor (40mm)
- 6 8-channel tip connector
- 7 Rack



Product portfolio

According to laboratory requirements, BioDolphin can be used with various types of vacuum pump / suction system like **Rocker 300 - BioDolphin**, **Protatable Suction System BioVac 225 - BioDolphin**, **Lafil 400 - BioDolphin**.




Rocker 300 - BioDolphin



Lafil 400 - BioDolphin



BioVac 225 - BioDolphin



199100-00	BioDolphin - LB Suction Kit (Light Blue)
199110-00	BioDolphin - DB Suction Kit (Deep Blue)
199120-00	BioDolphin - P Suction Kit (Pink)
199130-00	BioDolphin - G Suction Kit (Green)
197403-11	Lafil 400 - BioDolphin Protatable Suction System, AC110V, 60Hz
197403-22	Lafil 400 - BioDolphin Protatable Suction System, AC220V, 50Hz
167226-11	BioVac 225 - BioDolphin Protatable Suction System, AC110V, 60Hz
167226-22	BioVac 225 - BioDolphin Protatable Suction System, AC220V, 50Hz
167308-11	Rocker 300 - BioDolphin Vacuum Suction System, AC110V, 60Hz
167308-22	Rocker 300 - BioDolphin Vacuum Suction System, AC220V, 50Hz
199100-30	Handle
199100-46	1-channel tip adaptor
199100-47	1-channel S.S adaptor (40mm)
199100-48	1-channel S.S adaptor (80mm)
199100-49	8-channel S.S adaptor (40mm)
199100-50	8-channel tip connector
199100-61	Rack
197000-12	PES Suction Bottle(1200ml) with Spill- proof buoy
167300-11/22	Rocker 300 Vacuum pump, AC110V,60Hz / AC220V, 50Hz

ROCKER

6F., No.6, Zhulin Rd., Linkou Dist.,
New Taipei City 244, Taiwan (R.O.C.)
T/ +886-2-26033311
F/ +886-2-26036622
<http://www.rocker.com.tw>

Distributor