



## Supplementary Information of Rocker's Accessories

### Silicone Stopper Size



Part No.	Item	A Top	B Bottom	C Height	Hole
<b>167110-18</b>	#18 Silicone stopper	Φ 50.5 mm	Φ 44 mm	37 mm	Φ 9.2 mm
<b>167110-17</b>	#16 Silicone stopper	Φ 44 mm	Φ 37.5 mm	35.5 mm	Φ 9.2 mm
<b>167110-16</b>	# 8 Silicone Stopper	Φ 41 mm	Φ 33.4 mm	26 mm	Φ 9.2 mm
<b>167110-16-1</b>	# 8 Silicone Stopper	Φ 41 mm	Φ 33.4 mm	26 mm	Φ 15.6 mm
<b>167110-05</b>	# 5 Silicone Stopper	Φ 29 mm	Φ 24.5 mm	25.5 mm	Φ 9.6 mm





8號矽膠塞 (167110-16)  
No. 8 silicone stopper  
with 9.2 mm opening

16號矽膠塞 (167110-17)  
No. 16 silicone stopper  
with 9.2 mm opening

18號矽膠塞 (167110-18)  
No. 18 silicone stopper  
with 9.2 mm opening

Unit: mm

一般公差 TOLERANCE	±0.1/±1.0*	日期 DATE	比例 SCALE	料號 ITEM NO.
表面處理 SURFACE		材質 MATERIAL	視角 Protection	TEL
設計 DESIGNER	黃振榮	品名 DRAWING NAME	矽膠塞	FAX
審圖 INSPECT		部番 DRAWING NO.		07-2273358
				07-2273716
				洛科儀器股份有限公司 Rocker Scientific Co., Ltd.