

How to order?

You can select Combination A, B or C

A

Step 1. Choose the manifold




MultiVac 300 - MB
3 Places Stainless Steel Manifold
(180300-02)



MultiVac 600 - MB
6 Places Stainless Steel Manifold
(180600-02)

B

Step 1. Choose the manifold



MultiVac 300 - MS
3 Places Stainless Steel Manifold
(180300-01)



MultiVac 600 - MS
6 Places Stainless Steel Manifold
(180600-01)



Step 2. Choose the SF funnel



Rocker, 100 ml
SF Stainless Steel Funnel
(180100-01)



Rocker, 300 ml
SF Stainless Steel Funnel
(180100-03)



Rocker, 500 ml
SF Stainless Steel Funnel
(180100-05)



Stainless Steel Cover
for SF 100 ml funnel
(180100-11)



Stainless Steel Cover
for SF 300 ml / 500 ml funnel
(180100-13)



Step 2. Choose the filter holder



Rocker, 300 ml
MF 3, Magnetic Filter Holder
(200300-01)



Rocker, 500 ml
MF 5, Magnetic Filter Holder
(200500-01)



Rocker, 300 ml
VF 6, Glass Filter Holder
(167120-31)

to complete your order

C

Step 1. Choose the manifold



MultiVac 300
3 Places Stainless Steel Manifold Base
(180300-00)



MultiVac 600
6 Places Stainless Steel Manifold Base
(180600-00)



Step 2. Choose the disposable funnel and adaptor



Sartorius, 100 ml / 250 ml
Microsart® disposable plastic funnel
(16A07-10---N/
16A07-25---N)



Millipore, 100 ml / 250 ml
Microfil® disposable plastic funnel
(MIHAWG100 / MIHAWG250)



Pall, 100 ml / 250 ml
Sentino® disposable plastic funnel
(4870 / 4871)

+

+

+



SS316 Adaptor for Sartorius
disposable plastic funnel
(190100-74)



SS316 Adaptor for Millipore
disposable plastic funnel
(190100-73)



SS316 Adaptor for Pall
disposable plastic funnel
(190100-75)

* SS316 adaptor is only available for above-mentioned disposable funnels.

Customized Design

We can offer customized adaptor for your funnel

