

FR.RES.154 / B

Effective date: 05/01/2024

Disp&Flow - Celery (BIO.045.1)

General informations

Product	Rapid immunochromatographic test for qualitative determination of celery antigen in foods, kitchen and production facilities
Category number	BIO.045.1
Lot number	045.24.1
Mfg date	may-23
Expiry date	2025/05/31

Kit content

#	Item	Component	Qty	Units	Lot
1	BIO.045.1	Strips individually packaged into sealed pouch	1	pcs	DF240228-045
2	S005WZ	Specimen extraction buffer, 3ml	1	pcs	S005230904
3	IFU-BIO.045.1-FR	Instruction Biotem Test Celery, French	1	pcs	FRv20231004

Quality control procedure

Date	29/02/2024
Person	Brun Cindy
Result	Passed

Parameters of QC	Control material	Short description	Required range	Actual value
Sensitivity	QX392_904	Celery root extract	>=32	32
Max. dilution giving positive result				
Specificity	Negative Control	Vegetable extract (fennel, coriander, parsley)	Neg (1:100)	Neg

Product is compliant with specifications.

NOT INTENDED FOR HUMAN APPLICATION

Project Manager

Brun Cindy

Date : 29/02/2024



Quality Assurance

Gallois Marjorie

Date : 01/03/2024

