



15mL and 50mL PerformR® Centrifuge Tubes

PerformR® centrifuge tubes are designed for your everyday centrifugation and storage needs. These tubes include large white solvent resistant labeling areas and black easy-to-read graduations. They are available with either a flat cap that is easy to label or a plug style cap for horizontal applications. Always use appropriate cushion or adapters and biohazard containment precautions. Sterile packaging is sterilized with a validated and approved ISO 11137 method.

Racked tubes are packaged in a reusable, autoclavable (122°C) and freezable (-80°C) polypropylene rack that interlocks to allow you to create a workstation.

Polypropylene racks available separately and in colors, see page 120.



15mL PerformR® Centrifuge Tubes with Flat Caps

Item	Package	Qty/Pk	Sterile	Qty/Cs
3136-345	50/Bag, 10 Bags/Case	50	✓	500
3146-335	25/Plastic Rack, 12 Racks/Case	25	✓	300
3146-870	Bulk, Tubes & Caps in Separate Bags	500		500

15mL PerformR® Centrifuge Tubes with Plug-Style Caps

Item	Package	Qty/Pk	Sterile	Qty/Cs
3137-345	50/Bag, 10 Bags/Case	50	✓	500
3147-335	25/Plastic Rack, 12 Racks/Case	25	✓	300
3147-870	Bulk, Tubes & Caps in Separate Bags	500		500

Note: 15mL tubes don't work with Thermo (010-0378) centrifuge adapters



50mL PerformR® Centrifuge Tubes with Flat Caps

Item	Package	Qty/Pk	Sterile	Qty/Cs
3186-345	50/Bag, 10 Bags/Case	50	✓	500
3196-335	25/Plastic Rack, 12 Racks/Case	25	✓	300
3196-870	Bulk, Tubes & Caps in Separate Bags	500		500

50mL PerformR® Centrifuge Tubes with Plug-Style Caps

Item	Package	Qty/Pk	Sterile	Qty/Cs
3187-345	50/Bag, 10 Bags/Case	50	✓	500
3197-335	25/Plastic Rack, 12 Racks/Case	25	✓	300
3197-870	Bulk, Tubes & Caps in Separate Bags	500		500