



NOTES:

1. For reference only. This is not a quality or control document.
2. Wells are on a 9mm spacing/pitch in both directions. Well pattern is centered about SBS footprint.

THIS DOCUMENT CONTAINS INFORMATION WHICH IS THE SOLE PROPERTY OF LABCON N.A.. THIS DOCUMENT IS RECEIVED IN CONFIDENCE AND ITS CONTENTS MAY NOT BE DISCLOSED IN ANY FORM, BE THEY ELECTRONIC OR PAPER, WITHOUT THE PRIOR WRITTEN CONSENT OF LABCON NORTH AMERICA.

UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES.

INTERPRET DIMENSIONS AND TOLERANCE PER ANSI/ASME Y14.5M-1994

DOCUMENT SET IN THIRD ANGLE PROJECTION

DRAWN	DATE
V.Happ	10/27/2014

labcon North America

3700 Lakeville Hwy, Petaluma, CA 94954
tel (707) 766-2100 | fax (707) 766-2199 | www.labcon.com

TITLE Customer Service Drawing:
2.2mL Microplate

SIZE	DWG. NO.	REV.
A	CS-5495-1 (Microplate, 2.2mL)	0

SCALE 1:2	DO NOT SCALE DRAWING	SHEET 1 of 1
-----------	----------------------	--------------