



Metal Free® Centrifuge Tubes

For applications requiring the assurance
of low background metals



- Made of 100% medical grade ultra low metal virgin resins
- Certified below one part per billion (<1ppb) for trace metals
- Natural colorant-free flat caps

www.labcon.com

MetalFree® Special Features

15mL and 50mL Metal Free® Centrifuge Tubes

Pure, additive-free flat caps for easy labeling



Large, smudge-resistant writing panel for sample identification



Thicker tubes withstand higher g-forces and freeze to -80° without breakage




15mL TUBE

50mL TUBE



Max Freeze Indicator



Accurate, easy-to-read graduations & reverse scale



Caps & tubes feature recycling codes



Tubes and caps made with virgin, 100% medical grade CFR21 compliant ultra low metal resins

Labcon's Metal Free® centrifuge tubes are specifically processed to ensure that common trace metals and elements that can interfere with experiments are below an ICP-MS detection level of one part per billion (<1ppb).

Metals included in testing are: Calcium, Magnesium, Zinc, Iron, Copper, Manganese, Aluminum, Silicon, Nickel, Vanadium, Sodium, Phosphorous, Cobalt, Chromium, Potassium, Lithium, Selenium, Cadmium, Mercury, Arsenic, and Lead.



Freezable



Autoclavable with caps loosely attached

15mL Metal Free® Centrifuge Tubes

Item #	Package	Sterile	Qty/Pk	Qty/Cs
3134-345	Bag (50 tubes/bag, 10 bags/cs)	✓	500	500
3144-335	Rack (25 tubes/rack, 20 racks/cs)	✓	500	500

15mL Tube Dimensions: Length: 118mm | Inside diameter: 15mm | Outside diameter: 17mm

50mL Metal Free® Centrifuge Tubes

Item #	Package	Sterile	Qty/Pk	Qty/Cs
3184-345	Bag (50 tubes/bag, 10 bags/cs)	✓	500	500
3194-335	Rack (25 tubes/rack, 20 racks/cs)	✓	500	500

50mL Tube Dimensions: Length: 115mm | Inside diameter: 25mm | Outside diameter: 29mm

PACKAGING OPTIONS:



Recycled Cardboard Fiber Rack



Bag (tubes & caps attached)



Labcon North America
Petaluma, CA, U.S.A.

+1.707.766.2100
info@labcon.com
www.labcon.com



Metal Free® Centrifuge Tubes Quality & Compliance

Certificates are available to download by lot number at: www.labcon.com/certificates.html



Product details are current when published; subject to change without notice. Labcon, the Labcon logo, and variations on the foregoing are registered (®) or unregistered trademarks and service marks, of Labcon North America. © Copyright 2017 Labcon. All rights reserved.

REV 1_10.10.17