



Disposal Sets Overview

Sets for the safe disposal of liquid and explosive laboratory waste.

- ✓ Electrically conductive material to avoid explosion hazards.
- ✓ Sets with a practical lid or solid ball valve funnel and grounding cables.
- ✓ Safety canister with mechanical level control.
- ✓ Optionally with collecting tray to protect against leaks.
- ✓ All sets can be used immediately without additional accessories.

All set components can also be reordered individually:



Component	Funnel			Canister			Collection container	
	Type	Ball valve	Ball valve	Hinged lid	Safety canister	Safety canister	Safety canister	Collecting tray
Model	ARNOLD	ARNOLD	Marco	-	-	-	-	-
Feature	Lance (rem.)	with lance (220 mm)	without lance	mechanical level control	-	mechanical level control	with base insert	-
Content	-	-	-	10 litre	10 litre	20 litre	10 litre	-
Thread	S60/61	S60/61	S60/61	S60/61	S60/61	S60/61	-	-
Set	Part number	317 621	317 621	318 960	108 042	107 953	108 043	117 985
A	307 320	•			•			
B	307 315			•		•		•
C	307 416			•	•			•
D	307 321		•				•	