

Water for Injection

WFI-quality water by PAN-Biotech

Support your cell culture and biomanufacturing processes with PAN-Biotech Water for Injection (WFI) quality water. Our WFI-quality water meets highest industry standards (EP / USP) and is suitable for a broad variety of upstream and downstream applications.

PAN-Biotech`s Water For Injection (WFI) quality water is prepared by reversed osmosis, deionization, ultrafiltration and sterile filtration from a Water for Injection (WFI) system that **meets or exceeds current EP / USP specifications.**

Release Criteria

Sterility
Conductivity
pH
Endotoxin
TOC

Specifications

sterile
< 1.1 $\mu\text{S}/\text{cm}$
5 to 7
< 0.25 EU /ml
< 500 ppb

Release criteria for WFI-quality water are also based on EP / USP specifications and guarantee an ease in **regulatory compliance.**

PAN-Biotech WFI-quality water is not intended for parenteral administration or the preparation of sterile dosage forms of prescription compounding.

Contact us at
www.pan-biotech.de

PAN-Biotech GmbH | Am Gewerbepark 6 | 94501 Aidenbach



WFI-quality water is offered in **flexible packaging solutions** precisely meeting your demands.

3D Bags

Volumes from 50 liter to 1,000 liter

Bag material

Multilayer film
with EVOH barrier,
PE contact layer

Ports

Flexible number of ports
Flexible port position

Tubing Material

Silicone
C-Flex
TPE



Description

Water for Injection (WFI) Quality Water
Water for Injection (WFI) Quality Water

Size

200 L Bag
500 L Bag

Product Number

CT-99200B
CT-99500B

2D Bags

Volumes from 500 ml to 50 liter

Bag material

Multilayer film
with EVOH barrier,
PE contact layer

Tubing Material

Silicone
C-Flex

Bag tubing: PE

Description

Water for Injection (WFI) Quality Water
Water for Injection (WFI) Quality Water

Size

10 L Bag
20 L Bag

Product Number

CT-99010B
CT-99020B



PET bottles

Volumes - 500 ml and 1000 ml

Description

Water for Injection (WFI) Quality Water
Water for Injection (WFI) Quality Water

Size

500 ml
1000 ml

Product Number

CT-991500
CT-991000



Contact us at
www.pan-biotech.de

PAN-Biotech GmbH | Am Gewerbepark 6 | 94501 Aidenbach

