

## Amersham Shrimp Alkaline Phosphatase



Alkaline Phosphatase catalyzes the removal of 5'- and 3'-Phosphate groups from DNA, RNA, Nucleotides and Proteins.

Alkaline phosphatase is suitable for dephosphorylation of 5' and 3'-overhangs, recessed ends, blunt ends, and single nucleotides. Following incubation, alkaline phosphatase can be easily inactivated by heating.

- Bacterial alkaline phosphatase recombinantly expressed in *E. coli*.
- Active in PCR and restriction endonuclease buffers.
- Inactivated by heating at 75°C for 5 minutes.

Use together with Exonuclease ([/shop/molecular-biology/pcr-and-amplification/gel-extraction-and-clean-up/illustra-exonuclease-i-p-03102](#)) I for cleaning up PCR reactions prior to downstream application such as DNA sequencing.

### Unit Definition

One unit is the amount of the enzyme required to dephosphorylate 1 µg of linearized pUC57 DNA termini in 10 minutes at 37°C.

## Product Specifications ^

---

### SELECT YOUR MODEL

---

Amersham Shrimp Alkaline Phosphatase



Parameter	Amersham Shrimp Alkaline Phosphatase
Includes	10x reaction buffer: 200 mM Tris-HCl, pH 8.0, 100 mM MgCl <sub>2</sub> ; Dilution buffer: 50 mM Tris-HCl, pH 8.0
Activity	1 unit/ $\mu$ l
Purity	The enzyme is purified to apparent homogeneity and is free of all contaminating endonucleases, exonucleases and ribonucleases
Pack size	1000 units