

**Product code : L0040**

**Product name : Colcemid 10 µg/ml in PBS (Demecolcin)**

CAS Number	Components	Quantity in g/l
10035-04-8	Calcium Chloride Dihydrate	0.13300000
7791-18-6	Magnesium Chloride Hexahydrate	0.10000000
7447-40-7	Potassium Chloride	0.20000000
7778-77-0	Potassium Phosphate Monobasic Anhydrous	0.20000000
7647-14-5	Sodium Chloride	8.00000000
7558-79-4	Sodium Phosphate Dibasic Anhydrous	1.15000000
477-30-5	Demecolcine	0.01000000
WATER		990.20700000