

Product code : L0601

Product name : Earle's Balanced Salts Solution w/o Calcium w/o Magnesium

CAS Number	Components	Quantity in g/l
7447-40-7	Potassium Chloride	0.40000000
7647-14-5	Sodium Chloride	6.80000000
7558-80-7	Sodium Phosphate Monobasic Anhydrous	0.12200000
50-99-7	D-Glucose Anhydrous	1.00000000
34487-61-1	Phenol Red Sodium Salt	0.01100000
144-55-8	Sodium Bicarbonate	2.20000000
WATER		989.46700000