

Product code : L0605

Product name : Hanks' Balanced Salts Solution w/o Ca w/o Mg w/o Na Bicarbonate w/o Phenol Red

CAS Number	Components	Quantity in g/l
7447-40-7	Potassium Chloride	0.40000000
7778-77-0	Potassium Phosphate Monobasic Anhydrous	0.06000000
7647-14-5	Sodium Chloride	8.00000000
7558-79-4	Sodium Phosphate Dibasic Anhydrous	0.04788000
50-99-7	D-Glucose Anhydrous	1.00000000
WATER		990.49212000