

Salt Solutions Formulation

	L0605 Liquid mg/l	L0606 Liquid mg/l	L0607 Liquid mg/l	L0611 Liquid mg/l	L0612 Liquid mg/l	X0507 Liquid mg/l	X0509 Liquid mg/l	X0510 Liquid mg/l	X0513 Liquid mg/l	P0153 Powder mg/l	P0154 Powder mg/l
Hank's Balanced Salts Solution HBSS											
Calcium Chloride Dihydrate	/	185	/	/	185	/	1850	/	/	/	185
Magnesium Sulfate Anhydrous	/	97.67	/	/	97.67	/	976.7	/	/	/	97.68
Potassium Chloride	400	400	400	400	400	4000	4000	4000	4000	400	400
Potassium Phosphate Monobasic Anhydrous	60	60	60	60	60	600	600	600	600	60	60
Sodium Bicarbonate	/	350	350	350	350	/	/	3500	/	/	/
Sodium Chloride	8000	8000	8000	8000	8000	80000	80000	80000	80000	8000	8000
Sodium Phosphate Dibasic Anhydrous	47.88	47.88	47.88	47.88	47.88	478.8	478.8	478.8	478.8	47.88	47.88
D.C.*											
D-Glucose Anhydrous	1000	1000	1000	1000	1000	10000	10000	10000	10000	1000	1000
Phenol Red Sodium Salt	/	11	/	11	/	/	110	/	110	/	11

* Other Components