

Product code : L0611

Product name : Hanks' Balanced Salts w/o Calcium w/o Magnesium w/ Na Bicarbonate w/ Phenol Red

CAS Number	Components	Quantity in g/l
7447-40-7	Potassium Chloride	0.40000000
7778-77-0	Potassium Phosphate Monobasic Anhydrous	0.06000000
7647-14-5	Sodium Chloride	8.00000000
7558-79-4	Sodium Phosphate Dibasic Anhydrous	0.04788000
50-99-7	D-Glucose Anhydrous	1.00000000
34487-61-1	Phenol Red Sodium Salt	0.01100000
144-55-8	Sodium Bicarbonate	0.35000000
WATER		990.13112000