

**SUBJECT** : Adjustment of the pH specifications of L0615 : Dulbecco's Phosphate Buffered Saline w/o Calcium w/o Magnesium

**1- INFORMATION ON PRODUCT CONCERNED BY THE NOTIFICATION :**

Catalogue Number : L0615  
Batch Number : Not applicable  
Production date : Batches produced from 01/08/2022  
Request Number : Not applicable

**2- DESCRIPTION OF THE NOTIFICATION :**

Description : After analyzing more than 200 batches of L0615 produced, Biowest came to the conclusion that its pH specifications (7.0-7.6) were a little too low.  
The average pH of those batches is around 7.5. So the new specifications will be 7.2-7.8.

Documents linked : Technical Data Sheet and Safety Data Sheet of L0615

**3- PREVENTIVE / CORRECTIVE ACTION TO REALISE :**

Description : Not applicable  
This change has no impact on product quality as no change is made on its composition.

Written by : M. Varlet

Position : Quality Officer

Date : 28/07/2022

Signature :



**BIO WEST**  
SAS au capital de 884 400 Euros  
Rue de la Caille - BP 1  
49340 NUAILLE - France  
Tél : +33 (0)2 41 46 42 42  
Fax : +33 (0)2 41 46 42 50  
RCS Angers B 311 511 111  
N° Siret 311 511 111 111  
N° TVA intracommunautaire FR 2131151111

For more information, please contact :

BioWest  
Rue de la Caille  
BP11  
49340 NUAILLE

Tél : 02 41 46 42 42 Fax : +33 (0)2 41 46 42 50  
Fax : 02 41 46 38 90