
Salt Solutions Formulation

	Dulbecco's Phosphate Buffered Saline DPBS	L0615 Liquid mg/l	L0616 Liquid mg/l	X0515 Liquid mg/l	X0520 Liquid mg/l	P0750 Powder mg/l
Inorganic Salts	Calcium Chloride Dihydrate	/	/	/	1330	/
	Magnesium Chloride Hexahydrate	/	/	/	1000	/
	Potassium Chloride	200	/	2000	2000	200
	Potassium Phosphate Monobasic Anhydrous	200	144	2000	2000	200
	Sodium Chloride	8000	9000	80000	80000	8000
	Sodium Phosphate Dibasic Anhydrous	1150	795	11500	11500	1150

* Other Components