

Product code : L0630
Product name : Versene

CAS Number	Components	Quantity in g/l
6381-92-6	EDTA Disodium Salt Dihydrate	0.29200000
7447-40-7	Potassium Chloride	0.19300000
7778-77-0	Potassium Phosphate Monobasic Anhydrous	0.19000000
7647-14-5	Sodium Chloride	7.99500000
7558-79-4	Sodium Phosphate Dibasic Anhydrous	1.15000000
14431-43-7	D-Glucose Monohydrate	0.19800000
WATER		989.98200000