

Product code : L0909

Product name : Trypsin 0.25% in PBS w/o Calcium w/o Magnesium w/ Phenol Red

CAS Number	Components	Quantity in g/l
7447-40-7	Potassium Chloride	0.20000000
7778-77-0	Potassium Phosphate Monobasic Anhydrous	0.20000000
7647-14-5	Sodium Chloride	8.00000000
7558-79-4	Sodium Phosphate Dibasic Anhydrous	1.15000000
50-99-7	D-Glucose Anhydrous	1.00000000
34487-61-1	Phenol Red Sodium Salt	0.00200000
9002-07-7	Trypsin 1:250 (Non-ACF)	2.50000000
WATER		986.94800000