

Product code : L0930

Product name : Trypsin - EDTA 1X in solution w/o Calcium w/o Magnesium w/ Phenol Red

CAS Number	Components	Quantity in g/l
6381-92-6	EDTA Disodium Salt Dihydrate	0.22140000
7447-40-7	Potassium Chloride	0.40000000
7647-14-5	Sodium Chloride	8.00000000
50-99-7	D-Glucose Anhydrous	1.00000000
34487-61-1	Phenol Red Sodium Salt	0.00200000
9002-07-7	Trypsin 1:250 (Non-ACF)	0.50000000
144-55-8	Sodium Bicarbonate	0.58000000
WATER		989.29660000