

Product code : L0932

**Product name : Trypsin 0.25% - EDTA 0.02% in HBSS w/o Calcium w/o Magnesium w/
Phenol Red**

CAS Number	Components	Quantity in g/l
6381-92-6	EDTA Disodium Salt Dihydrate	0.20000000
7447-40-7	Potassium Chloride	0.40000000
7778-77-0	Potassium Phosphate Monobasic Anhydrous	0.06000000
7647-14-5	Sodium Chloride	8.00000000
7558-79-4	Sodium Phosphate Dibasic Anhydrous	0.04800000
50-99-7	D-Glucose Anhydrous	1.00000000
34487-61-1	Phenol Red Sodium Salt	0.01000000
9002-07-7	Trypsin 1:250 (Non-ACF)	2.50000000
144-55-8	Sodium Bicarbonate	0.35000000
WATER		987.43200000