

**Product code : L0940**

**Product name : Trypsin - EDTA 1X in PBS w/o Calcium w/o Magnesium w/o Phenol Red**

CAS Number	Components	Quantity in g/l
6381-92-6	EDTA Disodium Salt Dihydrate	0.25480000
7447-40-7	Potassium Chloride	0.20000000
7778-77-0	Potassium Phosphate Monobasic Anhydrous	0.20000000
7647-14-5	Sodium Chloride	7.95000000
7558-79-4	Sodium Phosphate Dibasic Anhydrous	1.15000000
50-99-7	D-Glucose Anhydrous	1.00000000
9002-07-7	Trypsin 1:250	0.50000000
WATER		988.74520000