

**Product code : P0120**

**Product name : Glasgow MEM BHK21 w/ L-Glutamine w/o Na Bicarbonate w/o Tryptose Phosphate Broth**

CAS Number	Components	Quantity in g/l
10043-52-4	Calcium Chloride Anhydrous	0.20000000
7487-88-9	Magnesium Sulfate Anhydrous	0.09767000
7782-61-8	Ferric Nitrate Nonahydrate	0.00010000
7447-40-7	Potassium Chloride	0.40000000
7647-14-5	Sodium Chloride	6.40000000
13472-35-0	Sodium Phosphate Monobasic Dihydrate	0.12400000
50-99-7	D-Glucose Anhydrous	4.50000000
1119-34-2	L-Arginine Monohydrochloride	0.04200000
30925-07-6	L-Cystine Dihydrochloride	0.03129000
56-85-9	L-Glutamine	0.29200000
5934-29-2	L-Histidine Monohydrochloride Monohydrate	0.02100000
73-32-5	L-Isoleucine	0.05240000
61-90-5	L-Leucine	0.05240000
657-27-2	L-Lysine Monohydrochloride	0.07310000
63-68-3	L-Methionine	0.01500000
63-91-2	L-Phenylalanine	0.03300000
72-19-5	L-Threonine	0.04760000
73-22-3	L-Tryptophan	0.00800000
69847-45-6	L-Tyrosine Disodium Salt Dihydrate	0.05219000
72-18-4	L-Valine	0.04680000
67-48-1	Choline Chloride	0.00200000
137-08-6	D-Ca Pantothenate	0.00200000
59-30-3	Folic Acid	0.00200000
87-89-8	Myo-Inositol	0.00360000
98-92-0	Nicotinamide (Nicotinic acid amide)	0.00200000
65-22-5	Pyridoxal Hydrochloride	0.00200000
83-88-5	Riboflavin	0.00020000
67-03-8	Thiamine Hydrochloride	0.00200000
34487-61-1	Phenol Red Sodium Salt	0.01600000
WATER		987.47965000