

**Product code : P0154**

**Product name : Hanks' Balanced Salts w/ Calcium w/ Magnesium w/o Sodium Bicarbonate w/ Phenol Red**

CAS Number	Components	Quantity in g/l
10035-04-8	Calcium Chloride Dihydrate	0.18500000
7487-88-9	Magnesium Sulfate Anhydrous	0.09768000
7447-40-7	Potassium Chloride	0.40000000
7778-77-0	Potassium Phosphate Monobasic Anhydrous	0.06000000
7647-14-5	Sodium Chloride	8.00000000
7558-79-4	Sodium Phosphate Dibasic Anhydrous	0.04788000
50-99-7	D-Glucose Anhydrous	1.00000000
34487-61-1	Phenol Red Sodium Salt	0.01100000
WATER		990.19844000