

Product code : P0350

Product name : Leibovitz L 15 Medium w/ L-Glutamine

CAS Number	Components	Quantity in g/l
10035-04-8	Calcium Chloride Dihydrate	0.18500000
7786-30-3	Magnesium Chloride Anhydrous	0.09366000
7487-88-9	Magnesium Sulfate Anhydrous	0.09767000
7447-40-7	Potassium Chloride	0.40000000
7778-77-0	Potassium Phosphate Monobasic Anhydrous	0.06000000
7647-14-5	Sodium Chloride	8.00000000
7558-79-4	Sodium Phosphate Dibasic Anhydrous	0.19000000
56-41-7	L-Alanine	0.45000000
56-40-6	Glycine	0.20000000
74-79-3	L-Arginine Free Base	0.50000000
70-47-3	L-Asparagine Anhydrous	0.25000000
7048-04-6	L-Cysteine Monohydrochloride Monohydrate	0.15717600
56-85-9	L-Glutamine	0.30000000
71-00-1	L-Histidine	0.25000000
73-32-5	L-Isoleucine	0.25000000
61-90-5	L-Leucine	0.12500000
657-27-2	L-Lysine Monohydrochloride	0.07500000
63-68-3	L-Methionine	0.15000000
63-91-2	L-Phenylalanine	0.25000000
56-45-1	L-Serine	0.20000000
72-19-5	L-Threonine	0.60000000
73-22-3	L-Tryptophan	0.02000000
60-18-4	L-Tyrosine	0.30000000
72-18-4	L-Valine	0.20000000
67-48-1	Choline Chloride	0.00100000
137-08-6	D-Ca Pantothenate	0.00100000
146-14-5	Flavin Adenine Dinucleotide Disodium Salt	0.00010000
59-30-3	Folic Acid	0.00100000
87-89-8	Myo-Inositol	0.00200000
98-92-0	Nicotinamide (Nicotinic acid amide)	0.00100000
58-56-0	Pyridoxine Hydrochloride	0.00100000
67-03-8	Thiamine Hydrochloride	0.00100000
59-23-4	D-Galactose	0.90000000
34487-61-1	Phenol Red Sodium Salt	0.01100000
113-24-6	Sodium Pyruvate	0.55000000
WATER		985.22739400