

Product code : P0390

Product name : McCoy's 5A w/ L-Glutamine w/o Sodium Bicarbonate

| CAS Number | Components | Quantity in g/l |
|------------|---|-----------------|
| 10035-04-8 | Calcium Chloride Dihydrate | 0.13246000 |
| 7487-88-9 | Magnesium Sulfate Anhydrous | 0.09766000 |
| 7447-40-7 | Potassium Chloride | 0.40000000 |
| 7647-14-5 | Sodium Chloride | 6.46000000 |
| 7558-80-7 | Sodium Phosphate Monobasic Anhydrous | 0.50435000 |
| 50-99-7 | D-Glucose Anhydrous | 3.00000000 |
| 56-40-6 | Glycine | 0.00750000 |
| 56-41-7 | L-Alanine | 0.01390000 |
| 1119-34-2 | L-Arginine Monohydrochloride | 0.04210000 |
| 70-47-3 | L-Asparagine Anhydrous | 0.04500000 |
| 56-84-8 | L-Aspartic acid | 0.01997000 |
| 52-90-4 | L-Cysteine Free base | 0.03150000 |
| 56-86-0 | L-Glutamic Acid | 0.02210000 |
| 56-85-9 | L-Glutamine | 0.21920000 |
| 5934-29-2 | L-Histidine Monohydrochloride Monohydrate | 0.02096000 |
| 51-35-4 | L-Hydroxy-L-Proline | 0.01970000 |
| 73-32-5 | L-Isoleucine | 0.03936000 |
| 61-90-5 | L-Leucine | 0.03936000 |
| 657-27-2 | L-Lysine Monohydrochloride | 0.03650000 |
| 63-68-3 | L-Methionine | 0.01490000 |
| 63-91-2 | L-Phenylalanine | 0.01650000 |
| 147-85-3 | L-Proline | 0.01730000 |
| 56-45-1 | L-Serine | 0.02630000 |
| 72-19-5 | L-Threonine | 0.01790000 |
| 73-22-3 | L-Tryptophan | 0.00310000 |
| 69847-45-6 | L-Tyrosine Disodium | 0.02610000 |
| 72-18-4 | L-Valine | 0.01760000 |
| 50-81-7 | Ascorbic Acid | 0.00050000 |
| 67-48-1 | Choline Chloride | 0.00500000 |
| 58-85-5 | D-Biotin | 0.00020000 |
| 137-08-6 | D-Ca Pantothenate | 0.00020000 |
| 59-30-3 | Folic Acid | 0.01000000 |
| 87-89-8 | Myo-Inositol | 0.03600000 |
| 98-92-0 | Nicotinamide (Nicotinic acid amide) | 0.00050000 |
| 59-67-6 | Nicotinic Acid | 0.00050000 |
| 150-13-0 | P-Aminobenzoic Acid (PABA) | 0.00100000 |
| 65-22-5 | Pyridoxal Hydrochloride | 0.00050000 |
| 83-88-5 | Riboflavin | 0.00020000 |
| 67-03-8 | Thiamine Hydrochloride | 0.00020000 |
| 68-19-9 | Vitamine B12 | 0.00200000 |

Product code : P0390

Product name : McCoy's 5A w/ L-Glutamine w/o Sodium Bicarbonate

| CAS Number | Components | Quantity in g/l |
|------------|--------------------------|-----------------|
| 70-18-8 | L-Glutathione (Reduced) | 0.00050000 |
| 91079-38-8 | Bactopectone | 0.60000000 |
| 34487-61-1 | Phenol Red (Sodium Salt) | 0.01020000 |
| WATER | | 988.04118000 |

DOMINIQUE DUTSCHER SAS