

Product code : P0390

Product name : McCoy's 5A w/ L-Glutamine w/o Sodium Bicarbonate

CAS Number	Components	Quantity in g/l
10035-04-8	Calcium Chloride Dihydrate	0.13246000
7487-88-9	Magnesium Sulfate Anhydrous	0.09766000
7447-40-7	Potassium Chloride	0.40000000
7647-14-5	Sodium Chloride	6.46000000
7558-80-7	Sodium Phosphate Monobasic Anhydrous	0.50435000
50-99-7	D-Glucose Anhydrous	3.00000000
56-40-6	Glycine	0.00750000
56-41-7	L-Alanine	0.01390000
1119-34-2	L-Arginine Monohydrochloride	0.04210000
70-47-3	L-Asparagine Anhydrous	0.04500000
56-84-8	L-Aspartic acid	0.01997000
52-90-4	L-Cysteine Free base	0.03150000
56-86-0	L-Glutamic Acid	0.02210000
56-85-9	L-Glutamine	0.21920000
5934-29-2	L-Histidine Monohydrochloride Monohydrate	0.02096000
51-35-4	L-Hydroxy-L-Proline	0.01970000
73-32-5	L-Isoleucine	0.03936000
61-90-5	L-Leucine	0.03936000
657-27-2	L-Lysine Monohydrochloride	0.03650000
63-68-3	L-Methionine	0.01490000
63-91-2	L-Phenylalanine	0.01650000
147-85-3	L-Proline	0.01730000
56-45-1	L-Serine	0.02630000
72-19-5	L-Threonine	0.01790000
73-22-3	L-Tryptophan	0.00310000
69847-45-6	L-Tyrosine Disodium	0.02610000
72-18-4	L-Valine	0.01760000
50-81-7	Ascorbic Acid	0.00050000
67-48-1	Choline Chloride	0.00500000
58-85-5	D-Biotin	0.00020000
137-08-6	D-Ca Pantothenate	0.00020000
59-30-3	Folic Acid	0.01000000
87-89-8	Myo-Inositol	0.03600000
98-92-0	Nicotinamide (Nicotinic acid amide)	0.00050000
59-67-6	Nicotinic Acid	0.00050000
150-13-0	P-Aminobenzoic Acid (PABA)	0.00100000
65-22-5	Pyridoxal Hydrochloride	0.00050000
83-88-5	Riboflavin	0.00020000
67-03-8	Thiamine Hydrochloride	0.00020000
68-19-9	Vitamine B12	0.00200000

Product code : P0390

Product name : McCoy's 5A w/ L-Glutamine w/o Sodium Bicarbonate

CAS Number	Components	Quantity in g/l
70-18-8	L-Glutathione (Reduced)	0.00050000
91079-38-8	Bactopectone	0.60000000
34487-61-1	Phenol Red (Sodium Salt)	0.01020000
WATER		988.04118000

DOMINIQUE DUTSCHER SAS