

Composition

P04-09500

MEM Eagle with EBSS
w: stable Glutamine
w: 2.2 g/L NaHCO₃

| | Components | mg/L |
|------------------|--|----------|
| Inorganic Salts | Calcium chloride x 2H ₂ O | 264.92 |
| | Magnesium sulfate | 97.67 |
| | Potassium chloride | 400.00 |
| | Sodium chloride | 6,800.00 |
| | Sodium dihydrogen phosphate x H ₂ O | 140.00 |
| Other Components | D(+)-Glucose | 1,000.00 |
| | Phenol red | 10.00 |
| Amino Acids | L-Arginine x HCl | 126.00 |
| | L-Cystine | 24.00 |
| | L-Glutamine stable | 435.08 |
| | L-Histidine x HCl x H ₂ O | 42.00 |
| | L-Isoleucine | 52.00 |
| | L-Leucine | 52.00 |
| | L-Lysine x HCl | 72.50 |
| | L-Methionine | 15.00 |
| | L-Phenylalanine | 32.00 |
| | L-Threonine | 48.00 |
| | L-Tryptophan | 10.00 |
| | L-Tyrosine | 36.00 |
| | L-Valine | 46.00 |
| Vitamins | D-Calcium pantothenate | 1.00 |
| | Choline chloride | 1.00 |

| | | |
|--|--------------------|----------|
| | Folic acid | 1.00 |
| | myo-Inositol | 2.00 |
| | Nicotinamide | 1.00 |
| | Pyridoxal x HCl | 1.00 |
| | Riboflavin | 0.10 |
| | Thiamine x HCl | 1.00 |
| | NaHCO ₃ | 2,200.00 |

DOMINIQUE DUTSCHER SAS