

	Leibovitz´s L-15 Medium Components	P04-27055 mg/L
Inorganic Salts	Calcium chloride anhydrous	140,00
	Magnesium chloride x 6H ₂ O	200,00
	Magnesium sulfate dried	97,67
	Potassium chloride	400,00
	Potassium dihydrogen phosphate	60,00
	Sodium chloride	8000,00
	<i>di</i> -Sodium hydrogen phosphate	190,00
Other Components	D(+)-Galactose anhydrous	900,00
	HEPES	0
	Phenol red	10,00
	Sodium pyruvate	550,00
Aminoacids	L-Alanine	225,00
	L-Arginine Base	500,00
	L-Asparagine x H ₂ O	250,00
	L-Cysteine	120,00
	L-Glutamine	0
	Glycine	200,00
	L-Histidine Base	250,00
	L-Isoleucine	250,00
	L-Leucine	125,00
	L-Lysine x HCl	93,75
	L-Methionine	75,00
	L-Phenylalanine	125,00
	L-Serine	200,00
	L-Threonine	300,00
	L-Tryptophane	20,00
L-Tyrosine	300,00	
L-Valine	100,00	
Vitamins	D-Calcium pantothenate	1,00
	Choline chloride	1,00
	Folic acid	1,00
	myo-Inositol	2,00
	Nicotinamide	1,00
	Pyridoxol x HCl	1,00
	Riboflavine-5'-phosphate sodium salt x H ₂ O	0,1075
	Thiamine monophosphate chloride	1,00
	Weight (g / L)	13,69

Supplemented with:

- L-Glutamine 300 mg/l
- Stab. L-Glutamine 510 mg/l
- 25 mM Hepes 5958 mg/l (NaCl reduced to 7500 mg/l)