

Composition

P04-32505

HBSS

w/o: Phenol red

w: 0.35 g/L NaHCO₃

	Components	mg/L
Inorganic Salts	Potassium chloride	400.00
	Sodium chloride	8,000.00
	di-Sodium hydrogen phosphate dried	47.88
	Potassium dihydrogen phosphate	60.00
	Calcium chloride x 2H ₂ O	185.44
	Magnesium sulphate anhydrous	97.67
Other Components	D(+)-Glucose anhydrous	1,000.00
	Phenol red	0.00
	NaHCO ₃	350.00