

Composition

P04-33050P HBSS
w/o: Ca and Mg
w/o: NaHCO₃

	Components	mg/L
Inorganic Salts	Potassium chloride	400.00
	Sodium chloride	8,000.00
	di-Sodium hydrogen phosphate dried	47.88
	Potassium dihydrogen phosphate	60.00
	Calcium chloride x 2H ₂ O	0.00
	Magnesium sulphate anhydrous	0.00
Other Components	D(+)-Glucose anhydrous	1,000.00
	Phenol red	10.00