

## Composition

**P04-34500**

**HBSS**

**w/o: Ca and Mg**

**w/o: Phenol red**

**w: 0.35 g/L NaHCO<sub>3</sub>**

	Components	mg/L
Inorganic Salts	Potassium chloride	400.00
	Sodium chloride	8,000.00
	di-Sodium hydrogen phosphate dried	47.88
	Potassium dihydrogen phosphate	60.00
	Calcium chloride x 2H <sub>2</sub> O	0.00
	Magnesium sulphate anhydrous	0.00
Other Components	D(+)-Glucose anhydrous	1,000.00
	Phenol red	0.00
	NaHCO <sub>3</sub>	350.00