

## Compositions

DPBS (10x)	Components	P04-37500 mg/L
<b>Inorganisc Salts</b>	Potassium chloride	2000,00
	Potassium dihydrogen phosphate	2000,00
	Sodiumchloride	80000,00
	Calcium chloride x 2H <sub>2</sub> O	1330,00
	Magnesium chloride X 6 H <sub>2</sub> O	1000,00
	<i>di</i> -Natriumhydrogenphosphat wfr.	11500,00

DOMINIQUE DUBOISCHER SAS