

Composition**P04-39500****EBSS****w/o: Phenol red**

| | Components | mg/L |
|-----------------|--|----------|
| Inorganic Salts | Potassium chloride | 400.00 |
| | Sodium chloride | 6,800.00 |
| | Sodium dihydrogen phosphate x H ₂ O | 140.00 |
| | Calcium chloride anhydrous | 200.00 |
| | Magnesium sulfate x 7H ₂ O | 200.00 |
| Other Component | D(+)-Glucose anhydrous | 1,000.00 |