

Composition**P04-49505****HBSS (10x)**
w/o: Phenol red
w/o: NaHCO₃

	Components	mg/L
Inorganic Salts	Potassium chloride	4,000.00
	Sodium chloride	80,000.00
	di-Sodium hydrogen phosphate dried	478.80
	Potassium dihydrogen phosphate	600.00
	Calcium chloride x 2H ₂ O	1,854.40
	Magnesium sulphate anhydrous	976.70
Other Components	D(+)-Glucose anhydrous	10,000.00
	Phenol red	0.00