

**Composition****P04-50500****HBSS (10x)**  
**w/o: Ca and Mg**  
**w/o: NaHCO<sub>3</sub>**

	Components	mg/L
Inorganic Salts	Potassium chloride	4,000.00
	Sodium chloride	80,000.00
	di-Sodium hydrogen phosphate dried	478.80
	Potassium dihydrogen phosphate	600.00
	Calcium chloride x 2H <sub>2</sub> O	0.00
	Magnesium sulphate anhydrous	0.00
Other Components	D(+)-Glucose anhydrous	10,000.00
	Phenol red	100.00