

	ITS I (100 x)	ITS II (100 x)	ITS IV(100 x)
	P07-03100, P07-03110	P07-03200, P07-03210	P07-03400, P07-03410
Insulin	0,5 mg/ml	1,0 mg/ml	1,0 mg/ml
Transferrin	0,5 mg/ml	0,55 mg/ml	0,55 mg/ml
Sodiumselenite	0,5 µg/ml	0,5 µg/ml	0,5 µg/ml
EBSS	8,34 mg/ml	8,34 mg/ml	8,34 mg/ml
BSA	0	0	50 mg/ml
Linoleic acid water soluble	0	0	5 µl/ml
Selenic acid	0	0	0