

	ITS I (100 x)	ITS II (100 x)	ITS IV(100 x)
	P07-03100, P07-03110	P07-03200, P07-03210	P07-03400, P07-03410
<b>Insulin</b>	0,5 mg/ml	1,0 mg/ml	1,0 mg/ml
<b>Transferrin</b>	0,5 mg/ml	0,55 mg/ml	0,55 mg/ml
<b>Sodiumselenite</b>	0,5 µg/ml	0,5 µg/ml	0,5 µg/ml
<b>EBSS</b>	8,34 mg/ml	8,34 mg/ml	8,34 mg/ml
<b>BSA</b>	0	0	50 mg/ml
<b>Linoleic acid water soluble</b>	0	0	5 µl/ml
<b>Selenic acid</b>	0	0	0