

**Composition****P07-03410 ITS solution IV (100x)**

	Components	mg/L
Inorganic Salts	Calcium chloride anhydrous	200.00
	Magnesium sulfate x 7H <sub>2</sub> O	200.00
	Potassium chloride	400.00
	Sodium chloride	6,800.00
	Sodium dihydrogen phosphate x H <sub>2</sub> O	140.00
	Sodium selenite	0.500
Other Components	BSA	50,000.00
	D(+)-Glucose anhydrous	1,000.00
	Linoleic acid water soluble	5.00 (µl/ml)
	Insulin	1,000.00
	apo-Transferrin human	550.00