

**Product code : P0750**

**Product name : Dulbecco's Phosphate Buffered Saline w/o Calcium w/o Magnesium**

CAS Number	Components	Quantity in g/l
7447-40-7	Potassium Chloride	0.20000000
7778-77-0	Potassium Phosphate Monobasic Anhydrous	0.20000000
7647-14-5	Sodium Chloride	8.00000000
7558-79-4	Sodium Phosphate Dibasic Anhydrous	1.15000000
WATER		990.45000000

DOMINIQUE DUTSCHER CAS